



Competency Verification Record

Philips MRX HeartStart™ Defibrillator

Name _____ ID# _____ Unit/Dept _____ Date _____

Give this sheet to your manager / NEC to place in your on-unit competency file

COMPETENCY

Trainer: Document employee competence by signing when the employee demonstrates competence for all expectations within the job function.

SECTION ONE: USE OF MRX™ DEFIBRILLATOR

- A. Demonstrates AED mode
- B. Identifies information displayed in monitor mode
- C. Interacts with waves
- D. Interacts with alarms
- E. Monitors ECG
- F. Demonstrates Manual mode
- G. Demonstrates pacing functions
- H. Maintains the MRx

Competency Verified By _____

SECTION TWO: STAFF WITH PEDIATRIC RESPONSIBILITY

- A. Locates pediatric pads
- B. Uses barrel connector to attach pads for neonatal patients
- C. Adjusts MRx into pediatric mode
- D. Sets therapy knob to 1-10 joules and adjusts to desired joules

Competency Verified By _____ or Not Applicable _____

SECTION THREE: STAFF WITH INTERNAL PADDLE RESPONSIBILITY

- A. Locates therapy cable port
- B. Plugs in internal paddles
- C. Verbalizes that paddles will only deliver 50 joules regardless of setting

Competency Verified By _____ or Not Applicable _____

TRAINING ONLY

SECTION FOUR: STAFF USING THE FRx™ AED FOR PATIENT TRANSPORT

Complete on-line CBL for Philips AED; no competency assessment needed